

For Easy, Time saving, User-friendly dental training & education  
To meet practical demands, ADEM has been evolving day by day.

### ORT-KIT-3

## Orthodontic typodont for mini-implant practice

For practicing teeth activation by orthodontic appliances with mini-implant anchorage.  
Combination of regular soft wax and transparent hard resin form with anatomical rooted two colored metal teeth.

This ensures efficient and safe training with using mini-implant.



Partially made of transparent hard resin portion for embedding mini-implant.




Anteriors and 1st. bicuspid portion is soft wax form.

Case : Class One

Upper and lower 1st. Bicuspid are extracted

Embedding mini-implant in the middle of upper left 2nd. Bicuspid and 1st. molar in order to have traction of anteriors with mini-implant anchorage.



For further inquiry:



**TEL: +81-75-256-3015 FAX:+81-75-221-8032**

**<http://www.trimunt.co.jp/>**

-Your partner to build a better dental training and education-