

For Easy, Time saving, User-friendly dental training & education To meet practical demand, ADEM has been evolving day by day.

# Implant training models

Two types are available depend on practical demand



IMP-015-UL
Edentulous Model
without gingivae



IMP-016-UL Dentulous Model with gingivae

Residual teeth

7 3 1 3 4 3 4 7 4 3 1 3

2

### For practical trainings of following cases

- Drilling
- implant embedding
- · implant embedding for overdenture
- · sinus lifiting

- Drilling
- · Implant embedding
- Sinus lifiting practicerance with Schneiderian membrane and septum
- · Implant embedding for overdenture
- · Incision & suture
- GBR

## **Mountable with ADEM Simple Manikin**





Lower Jaw can be fixed easily by two stage positioner.

## Simple training tool for incision and suture practice



### Incision and suture practice

Made from 3 different layers presuming the gingivae, periosteum and alveolar bone.



For further inquiry:



ΓEL: +81-75-256-3015 FAX:+81-75-221-8032 http://www.trimunt.co.jp/